

## » QUICK REFERENCE CARD EVONE®

**! This Quick Reference Card does not replace the Instructions for Use of Evone® or Tritube® !**

### Introduction on Evone

Mechanical ventilator Evone is intended to be used in sedated / anesthetized patients (>40 kg ibw), who require FCV® or Jet ventilation.

### FCV® MODE:

FCV® ventilates a patient with a controlled and stable insufflation and exhalation flow (by suctioning), between a set minimum tracheal pressure ((P)EEP) and maximum tracheal pressure (Peak). FCV® aims to create linear increases and decreases of intratracheal pressure. FCV® is used in elective situations and requires a cuffed airway and tubes with ~2 – 10 mm ID.

Compatible with: Tritube, single / double lumen endotracheal tubes (size 5 /CH35 or higher).

### JET MODE:

Single and double lumen jet ventilation (60-150 bpm) is used in elective procedures or to liberate a patient from ventilation and requires an open airway.

Compatible with: Tritube (with deflated cuff), Jet catheter or rigid bronchoscope.

### Materials

- Evone Control Unit
- Evone Cartridge ①
- Evone Airway Adapter ②
- Humid-Vent Filter Pedi straight (HME Filter) ③
- Evone Breathing Tubing ④ OR Conventional Tube Adapter (CTA) ⑥
- Tritube, jet catheter or rigid bronchoscope ⑤ OR a conventional adult endotracheal tube (single lumen at least 5 mm ID or double-lumen at least CH35) ⑦
- Empty syringe (20mL) to check cuff
- Syringe with 2-5 mL saline and ~15 mL air to purge lumens
- Manometer
- Lubricant spray (e.g. silicone)

### Materials for alternative ventilation Tritube:

- Ventrain® and manometer

OR

### All tubes:

- Conventional tube/(laryngeal) mask (in parallel with Tritube)
- Conventional (balloon) ventilation equipment

### Assembly with Tritube




Fig. 1 Assembly of the Evone Breathing System and Tritube

### Assembly with conventional tubes




Fig. 2 Assembly of the Evone Breathing System and conventional adult endotracheal tube

### Installation and set up

- 1 Switch on Evone.
- 2 Perform Startup checks successfully.
- 3 Patient set up menu: select gender and fill out length. Accept default settings or start with last used.
- 4 Check and, if required, adapt alarm limits.

Note that default settings are:

- FIO<sub>2</sub> 50%
- Inspiratory Flow 12 L/min
- I:E ratio 1:1.0
- Peak 15 mbar
- EEP 5 mbar

### Intubation with tritube

- 1 Inflate cuff of Tritube - check for leakage - deflate and wrap cuff around Tritube.
- 2 Patient with increased risk on secretions: ask to clear the throat by coughing and swallowing any secretions.
- 3 Induce anesthesia (TIVA).
- 4 Visually assess larynx and remove secretions if present.
- 5 Take out stylet Tritube, spray with lubricant and put stylet back.
- 6 Bend Tritube in curve required for intubation.
- 7 Remove stylet after the tip has passed the vocal cords.
- 8 Advance Tritube while turning to facilitate insertion.
- 9 Pull back to the position aimed for to avoid tracheal contact with the tip.
- 10 Inflate the cuff to be sure Tritube's tip is free from the tracheal wall.
- 11 Flush both lumen with air by syringe.
- 12 Fixate Tritube.



### Ventilation with Tritube

- 1 Connect Tritube to Evone (ventilation lumen and pressure lumen).
- 2 Optional: start ventilation with the cuff deflated to allow deepening of anesthesia (Jet mode).  
*Note that the airway is open (risk on aspiration).*
- 3 Start ventilation with the cuff inflated (25-30 mbar) in FCV<sup>®</sup> mode when anesthesia is optimized. A triangular pressure curve appears on the screen (Fig. 3).



Fig. 3 FCV<sup>®</sup> mode active

- 4 If needed adapt ventilation settings:
  - FiO<sub>2</sub> as preferred
  - EEP as preferred
  - Peak to adjust Tidal Volume
  - Inspiratory Flow to adjust Minute Volume

### Ventilation with conventional tubes

- 1 Induce anesthesia (TIVA).
- 2 Intubate patient as usual with tube of choice.
- 3 Oxygenate patient as preferred to allow deepening of anesthesia.
- 4 Connect tube to CTA of Evone when anesthesia is optimized.
- 5 Start ventilation in FCV<sup>®</sup> mode. A triangular pressure curve appears on the screen (Fig. 3).
- 6 If needed adapt ventilation settings:
  - FiO<sub>2</sub> as preferred
  - EEP as preferred
  - Peak to adjust Tidal Volume
  - Inspiratory Flow to adjust Minute Volume

### Handling obstructions

- 1 Stop ventilation.
- 2 Fiercely flush the pressure lumen and/or ventilation lumen with 2-5 mL saline followed by ~15 mL air.
- 3 In case secretions are still present in ventilation lumen, remove secretions using a suction catheter.  
*Note that the airway needs to be open.*
- 4 **Optionally:** slightly turn Tritube to avoid any tracheal wall contact and inflate cuff.
- 5 Purge lumen again with 2 mL saline followed by air.
- 6 Re-start ventilation.

### Sedation and relaxation

Because of the small lumen (high resistance) of the breathing circuit, coughing may result in tube dislocation and **spontaneous breathing is not possible**. In case of light anesthesia (indicated by e.g. irregular pressure curves, coughing, BIS>60, TOF>90%):

#### Optimize anesthesia

- Optionally deflate cuff of Tritube temporarily, to reduce trachea stimuli and allow breathing/coughing. Optionally disconnect CTA temporarily to allow breathing / coughing

### Individual optimization of FCV<sup>®</sup> ventilation based on lung compliance (optional)

- 1 Find 'Best EEP'
  - Change both EEP and Peak stepwise by 1-2 mbar; keeping driving pressure (Peak – EEP) constant.
  - Monitor tidal volume (V<sub>T</sub>): Increased V<sub>T</sub> indicates increased respiratory system compliance (C<sub>RS</sub>) and improved ventilation.
  - Choose EEP setting resulting in highest V<sub>T</sub> for similar V<sub>T</sub> choose lowest EEP for circulatory reasons.

*Note: Do not change settings too rapidly. Adequate application of the following steps requires equilibration periods of at least 1 minutes.*

- 2 Find 'Best driving pressure'
  - Change Peak pressure stepwise by 1-2 mbar.
  - Monitor V<sub>T</sub>
    - Per mbar increase of driving pressure, V<sub>T</sub> is expected to increase with value of C<sub>RS</sub>;
    - If V<sub>T</sub> increases over-proportionally, C<sub>RS</sub> will increase -> improved ventilation;
    - If V<sub>T</sub> increases under-proportionally, C<sub>RS</sub> will decrease -> ventilation not further improved.
  - Chose Peak setting resulting in highest C<sub>RS</sub>.

*Note: This step might lead to the application of higher tidal volumes than generally advised by common guidelines.*

- 3 Find 'Best flow'
  - Adjust flow depending on measured etCO<sub>2</sub> and/or PaCO<sub>2</sub>.
  - To reduce etCO<sub>2</sub> and/or PaCO<sub>2</sub>: increase inspiratory flow
    - Results in higher frequency with same V<sub>T</sub> and higher minute volume.
  - To increase etCO<sub>2</sub> and/or PaCO<sub>2</sub>: decrease inspiratory flow
    - Results in lower frequency with same V<sub>T</sub> and lower minute volume.

### Wakening the patient

- 1 Set FiO<sub>2</sub> as preferred.

#### Tritube

Wake patient using one of the two ventilation options:

- 2 With inflated cuff (e.g. in case of aspiration risk) in FCV<sup>®</sup> mode.
- 3 Gently wake patient (no shaking).  
Deflate cuff and extubate when patient awakes.



- 2 With deflated cuff in Jet mode (*risk on aspiration*).
- 3 Open airway required.
- 4 Adapt settings if required (e.g. lower driving pressure with higher frequency may reduce tracheal stimuli).

#### Conventional tubes

Wake the patient:

- 2 Disconnect CTA from tube allow waking up using preferred method of oxygenation.

*Note that spontaneous breathing is not possible when the CTA is connected to the conventional adult endotracheal tube.*

### Ventinova Medical B.V.

**A** Meerenakkerplein 7 **T** +31 (0)40 751 60 20  
5652 BJ Eindhoven **E** info@ventinova.nl  
The Netherlands